



**Certified Reference Material of purity of  
Methyl hydroxybenzoate of food antiseptic**



**Producer: National Institute of Metrology**

**Address: No.18, Beisanhuan donglu Beijing Tel: +86-10-64524710  
Fax: +86-10-64524716 e-mail: crm-service@nim.ac.cn**



The certified reference material is mainly used for the component detection of relevant antiseptic in the field of foodstuff, medicine etc. It is also used for the quality control in the medication manufacturing department and detection laboratory, the calibration of analytical instrument, the validation of analytical method etc.

国家标准物质资源共享平台  
www.ncrm.org.cn

## 1、Preparation

The certified reference material was prepared from the relevant Methyl hydroxybenzoate pure substance. According its physical and chemical properties, the certified reference material was gained by extraction, crystallization method and then purified. The certified reference material was dried and disposed after the qualitative and primary analysis which meets the demand.

## 2、Traceability and certification

The CRM purity is certified by freezing point depressing method. The traceability is ensured by using measuring method and measuring apparatus for metrology speciality.

## 3、Certified value and uncertainty

Sample number	Name	Certified value (%)	Relative expanded uncertainty(%) $(k=2)$
GBW(E) 100074	Methyl hydroxybenzoate	99.6	0.2

The uncertainty induced by the certified result of principle purity is in consideration in evaluation of uncertainty of this CRM.

## 4、Test of homogeneity and stability

According to the requirement of national criterion for primary certified reference materials, tests for this certified reference material were carried out through random sampling by LC-RAD method. The minimum sampling amount is 5 mg. F-test method was used for homogeneities test, No statistically significant differences among bottles were observed at this sample size.

Random sampling at regular intervals was made to test the stability of the certified reference material. The effective duration for use of the certified reference material is 1 year. The stability of the CRM are kept on motoring by maker during this period of this time, any change of the value will be informed to the customer in time.

## 5、Instructions for packing, storage and use

1. packing: the certified reference material is packed in brown screw-capped amber glass bottles each containing less than 200mg.

2. Storage and use: The certified reference material should be stored in cooling

国家标准物质资源共享平台  
www.ncrm.org.cn

国家标准物质资源共享平台  
www.ncrm.org.cn

temperature, lightproof , drying conditions. It should be dried in  $P_2O_5$  airer over 48 hours to eliminate moisture before use.

3. The CRM is poisonous and noxious substance, so it should be avoided inbreathing and contacting with skin. It should be airproof as soon as possible after unsealing, and should be stored in stated conditions.



## STATEMENT

1. The CRM is only limited in use of scientific research and analytical measurement. Any loss caused by improper use and storage by customer will not be respond by maker.
2. Please check the kind, number and package as soon as the sample arriving. The compensate do not relate to any loss except for the CRM.
3. The maker only answer for the intact certificate with CRM cachet of nim. Please keep the certificate appropriately.
4. Please contact with technical consultant section, if more information related to the use of the CRM is needed.

---

National Institution of metrology P.R. China

Address: No. 18 East Road North San Huan, Beijing Post Code: 100013

Tel for distribution: +86(10)64524710. Tel for technical consultation: +86(10)64524776

Fax: +86(10) 64524716. Website: [www.nim.ac.cn](http://www.nim.ac.cn) ; [www.ncrm.org.cn](http://www.ncrm.org.cn)

